

AUTOCLAVE

TECHNICAL SERIES-2

Steam Sterilization Validation



Sanjeev patel
sanjeev@gampservices.com

A Brief Overview of Steam Sterilization:

Sterilization: Validated process used to render a product free of living microorganisms including bacterial endospores.

Depyrogenation: Removal or inactivation of bacterial endotoxin.

Various types of steam sterilizer:

- ❖ Saturated Steam
- ❖ Water Immersion
- ❖ Water Cascade System
- ❖ Air-Steam Mixtures
- ❖ Gravity Air Displacement (unpacked materials sterilization)
- ❖ Vacuum air Displacement (Packed materials)

The basic principles for validation are:

1. Must use BIs to demonstrate lethality
2. Must use thermometric/ thermocouples
3. Cycle development and description of load patterns are pre-requisites
4. Can do time/temperature or F_0 approach for control
5. Calibrate thermocouples both pre and again post
6. Must include “worst case” conditions
 - Maximum and minimum loads/ patterns
 - One run of reduced cycle time / temperature
 - Cold start for at least one of three runs per load pattern

Microbiological aspects of Steam Sterilization:

-D Value ,Z Value , $T_b = 121$ Degree C & Lethality (F_0)

▪ What is the D value?

-Decimal reduction time - The time required at a certain temperature to kill 90% (eg..reduce population by log 1) of the organisms being studied. Thus, after an organism is reduced by 1 D, only 10% of the original organisms remain. Dependant on microbe and initial numbers. The D value of 1.5 means it takes 1.5minutes to reduce 1 log (to 10%) @121 Degree. D value of 2.0 means more resistant while a D value of 1 min means less resistant.

▪ What is a Z value?

-Temperature change required to produce a 1 log reduction in D value.

■ What is F_0 ?

-The number of minutes to kill a specified number of microbes with a Z value of 10°C at a temp of 121.1°C .

$$F_0 = \Delta t \sum 10^{\frac{T-121}{z}}$$

⊕

where Δt = time interval between measurement of T

T = temperature of the sterilized product at time t

z = temperature coefficient, assumed to be equal to 10

□

If we assume a sterilization lasting 15 minutes, constantly at 121°C , we obtain:

$$F_0 = 15 \times 10^{\frac{121-121}{10}} = 15 \times 10^0 = 15 \times 1 = 15 \text{ min } \textit{utes}$$

indeed according to the definition of F_0

If we assume sterilization lasts 15 minutes, constantly at 111°C , we instead obtain:

$$F_0 = 15 \times 10^{\frac{111-121}{10}} = 15 \times 10^{\frac{-10}{10}} = 15 \times 10^{-1}$$

$$F_0 = \frac{15}{10} = 1.5 \text{ min } \textit{utes}$$

Therefore, a 15-minute sterilization at 111°C is equivalent, in terms of lethal effect, to 1.5 minutes at 121°C ; this can be easily expected if $z=10$.

■ Overkill

- Use many more microbes than would find on items typically autoclaved. Negates the need to test sample for bio load before running the cycle.
- Use a sterilisation time exceeding what is necessary to kill a large number of microbes. Negates the need to determine D value of microbe.
- Overkill is generally defined as a 12-log reduction in bio load